**Manuale accessori
Accessories Manual
Manual accesorios
Manuel accessoires**

**ORION 40
ORION 40DS ORION
60
ORION 80**

1 Special arms (VA, VATVCC, WA1, WA2)

1.1 Structure and swing arm installation


For the installation instructions regarding the structure and the swing arm, make reference to the "Installation manual" attached to the Product. In case of presence of the coaxial cable, join the connectors at the time of installation.

1.2 WA1 and WA2 monitor stand installation

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POSITION RINGS AND SEEGER FROM THE BOTTOM

 TECNO-GAZ	Accessory manual	MA_001	01/09/15	EN
		Rev.0	Pag. 5 di 24	

See drawing 188

Align and fit the monitor stand (2) on the pin of the swing arm (1).

Fit the washers (3) (4) following the instructions and in the following order:

- Brass spacer washer (3);
- Teflon spacer (4);
- Brass spacer washer (3).

Fasten the snap ring (5) to the pin using suitable pliers.

Fasten the closing cap (6) using the three screws (7).

Take the cables out of the hole as shown on the drawing.

1.3 WA1 and WA2 monitor installation and adjustment

See drawing 158

Installation. Install the monitor (1) on the stand and tighten by means of the 4 knobs (2) and relative washers (3).

The structure is supplied with power cables with earth and signal cables with BNC terminals. Join the cables (4) coming out of the structure to the specific connectors.

At the customer's request, the structure can also be supplied without cables. To fit the cables inside the structure, follow the instructions in paragraph 1.3 Cable positioning.

Adjustment. Loosen the 8 screws (5) and adjust the support so as to centre it on the hub as shown on the attached drawing, to prevent monitor tilting problems.

By loosening the knobs (2), the horizontal position of the monitor can be adjusted, while by loosening the two knobs (6) the swinging movement can be adjusted.

The arm is sold calibrated to support a monitor weighing about 8 kg. If, after monitor installation, the arm does not stably maintain the set position, the balance will have to be adjusted as indicated in paragraph 1.4.

1.4 WA1 and WA2 cable positioning (optional)

WARNING – Before fastening the central block of the structure to the bar, fit the bar cables by running these inside the pin

See drawing 159

In order to insert the monitor cable, the two screws (2) will have to be loosened and the cover (1) of the central hub will have to be removed to gain access to the compartment.

Remove the cable from the holes (3) of the pin and wrap it around this twice, leaving it loose and not tight. Afterwards insert the cable inside the horizontal arm.

Remove the screws (4) and the plastic cap (5) to facilitate the insertion of the cable inside the side pin.

Remove the two opposite outer covers (6) levering with a flat-head screwdriver. Remove the two opposite protections (8) inserting a pin or a screwdriver inside the holes (7) and press on both sides at the same time.

Insert and run the cable inside the C-shaped lever (9) along the entire length of the swing arm and insert it in the monitor stand pin. Finally, remove the cable from the monitor stand hole and connect the connectors. If necessary, remove the closing cap (10).

Reposition the removed covers (6,8,10).

1.5 Calibration of swing arm and horizontal arm stop (VA, VATVCC, WA1, WA2)

WARNING – always insert the stop screw. Failure to do so could damage and break the cable.

See drawing 164

If, after monitor installation, the arm is not stable in the set position, regulate compression.

Loosen the screws (2) and remove the cover (1).

By means of the adjustment pin (4) rotate the lever upwards and load the spring, if the arm tends to drop, or downwards to unload the spring, if the arm tends to lift up.

The swinging movement of the arm can also be adjusted upwards.

The Product is delivered without any swinging limits. To reduce upward swinging, fit the pin in the ring nut (3) and rotate it downwards until it is in

the required position.

The swinging movement downwards cannot be changed.

After making adjustments, return the covering (1) to its original position.

The possibility also exists of selecting the stop position of the horizontal arm. To do this, select the stop position corresponding to one of the 5 slot screws (5) and remove the screw (5). Remove the hexagon socket head retention screw (6) and insert it in the position of the previously-removed screw. Insert the slot screw (5) in the position of the removed retention screw. It is possible to choose between 6 different stop positions.

1.6 Monitor and maximum weight

VA arm


The VA arm is designed to sustain a 15" to 24" monitor, with a weight that can vary between 2 and 10 kg

WA1 arm


The WA1 arm is designed to sustain a 15" to 32" monitor, with a weight that can vary between 4 kg and 20 kg

WA2 arm



The WA2 arm is designed to sustain two 15" to 24" monitors, with a weight varying between 2 and 16 kg altogether



1.7 Movement


VA ARM

VATVCC ARM

WA1 E WA2 ARMS

2 Cameras

2.1 Description of appliance

Versions

The Camera is available for the following models:

ORION 40/40DS, ORION 60, ORION 80. The system consists of removable cameras fitted to the head of the Product, remote control to manage the functions and, on request, LCD monitor to display the images. This system offers the possibility of choosing between two types of camera:

- SD camera
- FULL HD camera


The camera body is of the Plug & Play type, as long as the Product is ready to accept this.

The SD camera body can be fitted on all Products ready to accept the camera. The FULL HD camera can only be fitted on Products ready to accept the FULL HD camera.


2.2 Description of keyboard operation

Control panel


The camera operation control panel is available in versions with wall panel and box to be fitted to shelf. Camera functions are: on/off I/O key (1), image rotation (2), Zoom adjustment (3), focus adjustment (automatic and manual) (4), exposure adjustment (automatic and manual) (5), freeze image (6), white balance (7)

BOX



PANNELLO A PARETE



2.3 CAMERA installation and grip

Camera body
installation

Align the camera with the support flange, in accordance with the shape of the connector. Press the lock lever, fit the camera up tight (1) and release the lever (2).


To remove the camera body from the Product, press the lock lever with your finger (3) and remove the camera (4) downwards, Release the lever.


	CAUTION – Make sure the grip and camera are firmly fitted. Risk of falling.
---	--

Grip installation

To install the sterilisable hand-piece, align it with the camera body and push it upwards (5) and turn it until the 3 stops are fastened in the respective holes on the hand-piece (6).

To remove the hand-piece from the lamp, simply press the three stops (7) at the same time and pull the hand-piece (8) downwards.

2.4 CAMERA power connection

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The wall camera control board must be powered at 24VDC, connecting the connector terminals (+/-V, Earth) to the terminals of the power supply board.

The wall control must be connected to the control board inside the dome by means of the “OW” cable. Connect the “OW” cable of the Product to the “OW” connector of the wall control board.

Connessione video LCD

The LCD display screen must be connected to the Video/Audio cables of the Product.

Depending on the type of camera, the cables to be connected are:

 TECNO-GAZ	Accessory manual	MA_001	01/09/15	EN
		Rev.0	Pag. 12 di 24	

SD Camera: - 'COMPOSITE' video signal cable;
- audio signal cable (on request).

Full HD camera: - 'Y' (COMPONENT) video signal cable;
- 'Pb' (COMPONENT) video signal cable;
- 'Pr' (COMPONENT) video signal cable;
- audio signal cable (on request).

HDSDI camera: - 'COMPOSITE' video signal cable;
- audio signal cable (on request).

2.5 Technical specifications

FULL HD TVCC "SONY FCB-EV7100"	
Image device	1/2.8 type CMOS
Image effective elements	Approx. 2,38 Megapixel
Horizontal resolution	Full HD 1080/60p (1920 x 1080)
Lens	10x optical zoom, f=3.8mm (wide) to 38mm (tele), F1.8 to F3.4
Digital zoom	12x (120x with optical zoom)
Visual angle (H)	67° (wide) to 7.6° (tele)
Minimum distance objective	10mm (wide), 800mm (tele) (Default: 320mm)
Minimum illumination	1.4 lux
White balance	Auto, Indoor
AE control	Auto, Manual
Focusing	Full Auto, Manual
Video output	HD: Analog component Y/Pb/Pr / SD: VBS: 1.0 Vp-p (sync negative) Y/C
Operating voltage	6 to12 Vdc
Electric absorption	3.7 Watt

SD TVCC “SONY FCBEX20DP”	
Image device	1/3 type Super JAD II CCD
Image effective elements	Approx. 440.000 pixels
Lens	10x optical zoom, f=5.1mm (wide) to 51mm (tele), F1.8 to F2.1
Digital zoom	10x (120x with optical zoom)
Visual angle (H)	52° (wide) to 5.4° (tele)
Minimum distance objective	15mm (wide), 800mm (tele)
Minimum illumination	0.25 lux
White balance	Auto, ATW, Indoor, Outdoor, One-push, Manual
Focusing	Auto (sensitivity, normal, low), Manual
Video output	SD: Analog VBS: 1.0 Vp-p Y/C
Operating voltage	6 to12 Vdc
Electric absorption	6.5 Watt


3 Keyboard

3.1 Wall keyboard

Wall remote control panel

On request, the product can also be fitted with a wall control for managing the functions by means of a supplementary keyboard.

The functions of the wall keyboard are the same as those described in paragraph 5.2 of the operator’s manual.


3.2 Capacitive (touch) keyboard

Control panel on
Product

On request, a capacitive keyboard can be fitted on the Product, with touch technology. By touching with your finger on the surface of the keyboard, the following functions can be activated:

- ON and OFF I/O (1) with green indicator LED
- Adjustment of light intensity by dragging your finger over the bar or touching the sun symbol keys (2). The display of the level of set intensity is indicated by 5 blue LEDs.
- Selection of colour temperature 4500K and 5000K (3). The display of the setting is indicated by the lighting up of the corresponding blue LED. With the lamp off, the 4500K LED indicates the presence of power voltage in the Product.
- Start of the “Endoled” function letter E (4). The display of the set function is indicated by the lighting up of the corresponding green LED. This function can only be used when the lamp is off.
- Adjustment of the light range (5). The keys extend or reduce the lit diameter. Display is by means of the lighting up of one of the 2 blue LEDs.

4 Infra-red remote control

4.1 Standard remote control

Functioning

On request, an IR receiving system can be installed where the central Product hub is located.


By pointing the remote control in the direction of the IR receiver, the controls of the Product keyboard can be operated

Control panel

The remote control, has a keyboard that permits making the same adjustments as the Product keyboard.

By means of the L key, the dome to be adjusted can be selected.


Screen

The remote control features a 3.5" LCD screen.

On the right is the selected dome (2). The 5 levels (3) indicate the set light

intensity.


In the lower part are indicated: the colour temperature (4), the selected diameter (5) the state of the Endoled function (6).

Maintenance

In the event of the Product not responding to the controls, make the following checks:

- Make sure the remote control is pointing towards the IR receiver placed in the central hub of the structure;
- Make sure the batteries are charged. To replace them, open the rear door by loosening the screws and replace them with new type AA ones.


4.2 Universal remote control

Control keyboard

A universal remote control is available which controls the lamp and camera

functions. For correct operation point the remote control in the direction of the central hub of the structure.

The keys perform the following functions:

- (1) lamp and camera switch-on and switch-off;
- (2) +/- adjustment of the selected function (lamp or camera);
- (3)(4) selection of unit to be adjusted (dome/satellite/camera);
- (4) in camera mode, the relative functions can be adjusted (keys 8);
- (5) cyclic change of colour temperature;
- (6) cyclic change of the diameter of the lit field;
- (7) Endoled function only available with lamp off;
- (8) camera functions (rotation, zoom, focus freeze, autoexposure, white balance).

5 Battery systems

5.1 Batteries for single lamps

	CAUTION – before using the Product, connect the faston to the battery lead connection. There could be a slight electrical discharge on the faston.
	CAUTION – before use, set the fuse on the terminal board.
	CAUTION – When disconnecting the Product from the power supply, the switch must be in off position (0) to prevent discharging the batteries.
	CAUTION – to ensure correct battery operation, replace them at least every 3 years.
	CAUTION – in case of used up battery charge, the lamp could go off.

Storage precautions The lamps are supplied with a series of two charged batteries.
Store the batteries at a temperature between -20°C and +40°C.
During storage, recharge the batteries at least once every 6 months.
A battery also ages during storage, so it is always best to use it as quickly as possible.

Operating precautions It is best to perform a test every month to check the operation of the exchange relay and the battery charge.
To do this, disconnect the Product from the supply mains and make sure switchover to the battery occurs automatically (the lamp must not go off).
When the Product is not used for more than one month, disconnect the batteries.
The batteries are an integral part of the power circuit and the switch from mains power to battery occurs automatically in case of a power interruption. When the mains power supply returns, the switch from battery to mains is automatic. A battery charger keeps the batteries charged.

 TECNO-GAZ	Accessory manual	MA_001	01/09/15	EN
		Rev.0	Pag. 19 di 24	

5.1.1 Batteries for ORION40/40DS

Battery specifications

Nominal voltage 12 V;
Capacity 12 Ah;
Maximum charge current 3 A;
Weight 3.75 kg;
Dimensions 151x98x94 (mm);
Operating range -20°C / 40°C.

Battery life

The batteries have a life of about 4 hours.
The batteries are subject to a number of charge and discharge cycles.

- If the batteries are completely discharged the number of theoretical life cycles is 250.
- If the batteries are discharged by 50% and subsequently recharged, the number of theoretical life cycles is 550.
- If the batteries are discharged by 30% and subsequently recharged, the number of theoretical life cycles is 1200.

Indicator operation

The battery charge is indicated by a coloured LED.
If the light is green, this means the batteries are charged, if it is orange, the batteries are being charged and if it is red this means the batteries are down.
When the indicator is red, connect the lamp to the power mains.

5.1.2 Batteries for ORION60 ORION80

Battery specifications

Nominal voltage 12 V;
Capacity 12 Ah;
Maximum charge current 6.8 A;
Weight 8.5 kg;
Dimensions 166x175x125 (mm);
Operating range -20°C / 40°C.

Battery life

The batteries have a life of about 3 hours.
The batteries are subject to a number of charge and discharge cycles.

- If the batteries are completely discharged the number of theoretical life cycles is 200.
- If the batteries are discharged by 50% and subsequently recharged, the number of theoretical life cycles is 400.
- If the batteries are discharged by 30% and subsequently recharged, the number of theoretical life cycles is 1200.

The system does not feature battery charge indicators.

 TECNO-GAZ	Accessory manual	MA_001	01/09/15	EN
		Rev.0	Pag. 20 di 24	

5.2 Batteries for small double lamps

Validity	This system is used on the ORION40+40 / 40DS+40DS
Product	The system consists of a power supply unit made up of two separate boxes, each of which controls a lamp, placed in communication with each other by means of a common cable. This cable must be connected to the respective terminals before commissioning in accordance with the numbering and wiring diagram. Each box contains two batteries. The system does not feature battery charge indicators.
Battery specifications	Nominal voltage 12 V; Capacity 12 Ah; Maximum charge current 3 A; Weight 3.75 kg; Dimensions 151x98x94 (mm); Operating range -20°C / 40°C. The system does not feature battery charge indicators.
Battery life	The batteries have a life of about 4 hours. The batteries are subject to a number of charge and discharge cycles. - If the batteries are completely discharged the number of theoretical life cycles is 250. - If the batteries are discharged by 50% and subsequently recharged, the number of theoretical life cycles is 550. - If the batteries are discharged by 30% and subsequently recharged, the number of theoretical life cycles is 1200.

5.3 Battery unit in cabinet (GBA)

5.3.1 General information

Validity	This system is used on the ORION 60+40, ORION 60+60, ORION80+60, ORION 80+80
Product	The system consists in a power supply unit for the lamps with battery charge system. The GBA allows the lamps to operate both from the mains at 230V or, in case of a power breakage, from a battery at 24 V. The switch from mains to battery is automatic in case of sudden power break or if the cabinet door is opened.


To return to mains operation, in the first case the mains supply must be available or in the second case, the door of the cabinet must be closed again. In both cases, the GBA switches over automatically.

The system features a voltmeter and an ammeter. The voltmeter indicates the battery charge while the ammeter indicates the charge current of the battery charger.

	CAUTION – Battery operation is limited to battery life. Make sure the voltage indicated by the voltmeter does not fall below 19V, otherwise the light could go off.
---	--

5.3.2 Battery charger charging cycle


Quick charge

The charge current is maximum and at the same time the battery voltage is below the timer interruption level. The battery charger LED will be red/orange.

The ammeter will show a current of 5A.

Timer

The battery charger is in mode with time control. The current of the battery charger is below maximum. The battery is usually charged between 80% and 95% when the mode with time adjustment starts (the LED becomes yellow). The battery voltage coincides with the quick charge interruption level. The battery charger remains in this mode until the time interval has passed.


The ammeter shows a decreasing current from 5A to 0.8A approx.

Charge in standby mode

The battery charger is in standby mode. The battery is completely charged (the LED becomes green). The charge voltage is at standby level, which means the battery charger can continue to be connected to the battery. The battery charger can return to quick charge if the battery is down.

Charge graph

Battery charger capacity	Charge current	Charge voltage	Stability voltage	Timer	Protection	
20-200Ah	5 A	29.4 V	27.4 V	4h ± 30min	Fuse	230V
20-200Ah	2 A	29.4 V	27.4 V	9h ± 30min	Fuse	110V

5.3.3 Parameters and indicators

On the front of the GBA the following information can be displayed:

Parameters	Battery voltage displayed by means of voltmeter. Battery charger charge current by means of ammeter.
Indicators	Green indicator showing supply of mains power voltage. Green indicator showing main lamp (1) operating. Green indicator showing satellite lamp (2) operating. Battery charger indicator with colours indicated at para, 5.2.2.

5.3.4 Maintenance

CAUTION


The GBA does not require any particular routine maintenance.
Clean it inside once a year using compressed air to remove dust from the component parts and make sure the connections are tight.

5.3.5 First switch-on

CAUTION

The batteries are dispatched charged and disconnected with the battery fuse (F13) attached separately.
The first step is to connect the batteries by means of the faston.
Position the battery fuse.
Connect the GBA to the power mains and then the lamps to the GBA as shown on the wiring diagram.
Before energizing the system, make sure all the fuses are engaged.
Energize the system and make sure the GBA is working properly in both the charging and discharging cycle.

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5.3.7 Battery life

CAUTION – to ensure correct battery operation, replace them every 3 years.

LAMP MODEL	LIFE (h)
ORION 60+60	2
ORION 60+40	2.5
ORION 80+80	2
ORION 80+40	2.5

At the end of the indicated operating span, the batteries must be recharged. To completely recharge the batteries, 6 hours are needed.

The batteries have a limited number of charge and discharge cycles.

- If the batteries are completely discharged, the theoretical number of operating cycles is 200.
- If the batteries are discharged by 50% and then recharged, the number of theoretical operating cycles is 400.
- If the batteries are discharged by 30% and then recharged, the number of theoretical operating cycles is 1200.

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Copyright