MultiSteril FAST

ITALIANO	Manuale Tecnico
ENGLISH	Technical Manual
FRANCAIS	Manuel Technique
ESPAÑOL	Manual técnico

0ZVMT0001
16/01/2018

INDICE

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INDEX

01	GENERAL WARNINGS	2
02	GENERAL INFORMATION	3
	<i>TECHNICAL FEATURES</i> <i>DIMENSION</i>	
03	DESCRIPTION OF COMPONENTS	5
	<i>CONTROL PANEL</i>	
04	INSTALLATION	9
	<i>PLACE PREARRANGEMENT</i> <i>INSTALLING THE WATER SUPPLY</i>	
05	DIAGRAMS AND OPERATION	12
	<i>PNEUMATIC DIAGRAM</i> <i>WIRING DIAGRAM</i> <i>OPERATION PHASES</i> <i>TECHNICAL MENU'</i>	
06	ERRORS	27
	<i>ERR 005</i> <i>ERR 007</i> <i>ERR 008</i> <i>ERR 009</i> <i>ERR 010</i> <i>ERR 011</i>	
07	OPERATING PROCEDURES	32
	<i>P1 - OPENING OF THE TOOL FRONT PANEL</i> <i>P2 - OPENING OF THE PLASTIC COVER OF THE LID</i> <i>P3 - REMOVAL OF THE REAR PANEL</i> <i>P4 - REMOVAL OF THE POWER BOARD</i> <i>P5 - TUB REMOVAL</i> <i>P6 - CLEANING OF THE FILTER OF THE CHARGE SOLENOID VALVE</i> <i>P7 - REPLACEMENT OF THE DRAIN PUMP</i> <i>P8 - DISASSEMBLY OF THE CHARGE SV</i> <i>P9 - DISASSEMBLY OF THE SENSOR LEVEL</i> <i>P10 - DISASSEMBLY OF THE PERISTALTIC PUMP</i>	
08	TESTING PROCEDURE	41

01**GENERAL WARNINGS**

TECNO-GAZ S.p.A. is responsible for the safety, reliability and performance of the equipment, if:

- installation, any modifications, calibrations or repairs are carried out by authorized technical staff using original spare parts.
- the electrical installation of the related local corresponds to current rules.
- the device is used according to the instructions in the user manual DVMA034.

Further notes:

- Ensure that the unit is powered with a suitable electrical system with grounding and with right voltage marked in the silver label.
- Do not remove the silver label.
- Before any work unplug the power cable from the socket.
- Before installation check the integrity of the device in order to verify any transport damage.
- Use only original spare parts.


02

GENERAL INFORMATION

MULTISTERIL is a multifunctional appliance able to implement automatically without any operator's intervention some of the main phases of the sterilization cycle, such as disinfection / clearing / rinse/ drying.


TECHNICAL FEATURES	
MECHANICAL	
<i>Working temperature</i>	+5°C ÷ +30°C
<i>Humidity at 30°C</i>	80%
<i>Dimensions (LxHxP)</i>	(415 x 415 x 530) mm
<i>Height with lid open</i>	750 mm
<i>Dept with lid open</i>	530 mm
<i>Weight (with empty tank)</i>	20 kg
ELECTRICAL	
<i>Power voltage</i>	230 V a.c. +/-10 % single phase
<i>MAX absorbed power</i>	1 kW
<i>Frequency</i>	50 / 60 Hz
<i>Power cord</i>	2 + 1 x 1mm ²
<i>Fuses</i>	6,3 A (6 x 20)mm
ULTRASOUND	
<i>N° transductor</i>	4
<i>Frequency</i>	39 kHz
WATER SUPPLY	
<i>Water Pressure</i>	1.5 ÷ 2 Bar
<i>MAX Temperature</i>	+30°C
<i>Dimension loading pipe</i>	3/4 G - lenght 3 m
TUB	
<i>Tub Volume</i>	9 l (Inox AISI 316)
<i>Volume water loading</i>	6 l
DETERGENT TANK	
<i>Volume</i>	1.5 l + 0.5 l (riserva)
<i>Material</i>	Polietilene
FILTERS	
<i>Tub filter</i>	Inox AISI 304
<i>Cartridge filter</i>	Inox AISI 304 (19 x 40) mm

DIMENSION


03

DESCRIPTION OF COMPONENTS


File 7ZVMA0001

A	Drying fan	B	Safety thermostat
C	Piezoelectric transducers	D	Detergent tank
E	Drain pump	F	Detergent pump
G	Charge solenoid valve	H	Power socket with ON/OFF switch
I	Silicone resistances	J	Safety microswitch
K	Probe PT100	S1	Overflow detector
S2	Level sensor		


File 7ZVMA0002

L	Power board	M	Detergent tank drain
N	Water charge EV	O	Water drain EV
P	Drain tube	Q	Gas spring
R	Drain tube support		

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

Do not obstruct lateral and rear vent grids.
Leave a space of at least 5 cm around the
appliance to ensure a correct ventilation

Position the appliance far
from heat sources or
strong electromagnetic
radiations

Minimum difference in height between machine
drain and local drain = 20cm


Working pressure = $1.5 \text{ bar} < \Delta P < 2 \text{ bar}$

Power socket = SCHUCO

INSTALLING THE WATER SUPPLY


01

Take the machine out of the package by
means of the lifting belts
(this operation should be carried out by two
operators)

02

Position the appliance on a perfectly horizontal
surface, able to support a weight of at least 35 kg.


OK

NO


03

Take the accessories kit out of the tank inside

04


Unscrew the cap from the charge coupler and remove the drain coupler cap

05

Appliance drain system connection to the local net



WARNINGS = minimum external diameter of the drain coupler 22 mm

SIFONE
SIFONE
SIPHONS

06

Connection of water charge hose and WATERBLOCK device (supplied and compulsory).

Insert the two seals provided into the corresponding threaded couplers before screwing them.


Filettatura rubinetto = 3/4 G

OPTIONAL:
The WATER BLOCK
has to be mounted
vertically between the
threaded ring nut of the
charge hose and the
cock (check that this is
set on 4)


07

Connection of detergent tank drain and connection to the power socket


- i** The assembling shall be performed by qualified operators, who must follow the indication provided by the use instructions and the quick guide coming with the device
- i** The technician who performed the assembling must take care of the operators' training for the usage of the installed machinery
- i** *AT THE END OF THE DAY OR IF THE MACHINE IS NOT USED FOR SOME DAYS CLOSE THE CHARGE WATER COCK*

PNEUMATIC DIAGRAM

File 4ZVMP0002

WIRING DIAGRAM

File 4ZVME0002

V	VESSEL
C	COVER
UD	CLEANER NOZZLE
UA	WATER NOZZLE
SD	CLEANER TANK
IA	WATER INLET
ID	CLEANER INLET
SA	WATER EXHAUST
TM	KEYBOARD MEMBRANE
PP	PERISTALTIC PUMP
LV2	TANK LEVEL
SL1	OVERFLOW DETECTOR
SL2	SENSOR LEVEL
ST	TEMPERATURE SENSOR
RC	ARMoured RESISTANCE
TS1/TS2	SAFETY THERMOSTAT (RC)
R1/R2/R3	RESISTOR SILICONE R1/R2/R3
EV1	CHARGE SOLENOID VALVE
TP1/TP2/TP3/TP4	TRANSDUCERS
M1	MICROSWITCH
F1/F2	FAST FUSE F6A
S1	POWER PLUG

OPERATION PHASES

After installing the tool, open the cock connecting the charge hose and switch on the appliance rotating the switch on the socket **(H)** to position – I – .
Upon the switch ON the logo and the appliance firmware version are displayed.
After some seconds the cycle selection page is shown.

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
- ▶ CICLO A
- CICLO B
- CICLO C


With lid lifted and tool ON lift the cap **(T)** for detergent refilling.
Pour detergent (the use of a funnel is suggested) until an acoustic/visual alarm informs that the level probe inside the tank **(D)** reaches the MAX level.

WARNING: Do not fill the detergent tank with tool OFF since the level probe (which is not power supplied), will not be able to indicate the tank level.


Insert the tools in to the net basket and then this in to the tank, being careful not to create contact between the walls and the basket (any contact between tank and metal could lead in time to breakage of the tank due to the vibrations of the ultrasounds).

Select the cycle pressing **SELECT** to scroll the menu vertically

Press **START**


PHASE 1 : WATER LOADING

As first operation, the tool loads the tank (**W**) with tap water by opening the charge valve (**G**) to allow flow to the nozzles (**Z**) located on the lid.

Water is loaded until the probe (located inside the level gauge block (**Q**)) reaches the maximum level (6L) : the control board then closes the charge solenoid valve.

2.5 min
(con $P_{\text{ingresso}} = 2.5$ bar)

PHASE 2 : DETERGENT LOADING


The peristaltic pump loads the detergent through the injector, in the percentage selected by the preset cycle.


Percentage selected	Volume (ml)	Cycles ex.*
1%	60	≈25
2%	120	≈12
3%	180	≈8
4%	240	≈6
5%	300	≈5
6%	360	≈4

*With a full tank (1500 ml)


%	sec
1	80
2	160
3	240
4	320
5	400
6	480

PHASE 3 : HEATING / DEGAS

	The silicone resistances on the tub start heating the water volume at the temperature set in the current cycle, temperature is then measured by the temperature probe PT100. (If it doesn't reach the temperature, the heating continues also during the washing -	5 min¹
 LAMPEGGIANTE	The device activates intermittently ultrasounds. This process reduces the gases dissolved in the liquid by increasing the efficiency of the next washing.	5 min¹

The heating and degassing stages, if set both, are carried out simultaneously for a total duration of 5 minutes.

PHASE 4 : ULTRASONIC CLEANING


The continuous activation of transducers **(C)** located on the tub bottom carries out a cleaning action. In this phase both solenoid valves (charge and drain) are closed. This phase lasts for the selected time.

from 1 min
to 30 min

SWEEP FUNCTION

If selected, this function cyclically changes the working intervals of transducers thus optimizing the tub cleaning.


PHASE 5 : DRAIN

Once the ultrasound phase has been completed, start the drain pump to allow the water-detergent solution draining from the tub.
The drain phase is over after 1.5 minute

≈ 1.5 min

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
LOADED WATER - Filling the tank through water inlet nozzles

ULTRASONIC RINSE - Activation of the ultrasound for 1 minute

EXHAUST

≈ 5 min

PHASE 7 : DRYING

The last phase is drying. The fan **(A)**, generates an air flow passing through the armored resistance **(S)**. Hot and dry air touches the instruments and dries them.


The appliance sucks fresh air from the slots in the rear area of the cover panel and ejects hot air from side slots on the lid.

Do not obstruct the air suction and vent slots.

During the 5 last minutes of the drying phase the resistance turns off while the fan still operates to allow for the cooling of the resistance and of the instruments inside the tub.

At the end of the cycle the appliance sends out three acoustic signals and displays the end of cycle.

15 min

During the drying phase or when ventilation is ON, do not switch OFF the appliance. In case of power cut off for any reason, immediately open the lid to avoid any damage to the ventilation and to the appliance lid.

ATTENTION

- *Do not open the lid when the appliance is operating. In case of accidental opening, the appliance stops operating. Once the lid is closed the appliance resumes its cycle.*
- *During the cleaning or rinsing phases the opening of the lid might cause water overflow.*
- *During the rinsing phase the lid opening stops the heating phase, however the fan goes on operating to reduce the resistance temperature. In this case wait for some minutes before handling the basket and the instruments.*

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
TECHNICAL MENU

To enter the technical menu switch On the tool and after the page with the TECNO-GAZ logo is displayed wait for the next page indicating that the tool is operating. Hold down the STOP button and at the same time press the OK button twice to enter the technical menu.

SELECT	OK	DEVICES TUB : water level into the tub DETERGENT: detergent tank level PT100°C : temperature measured by the probe on the tank <i>From peripherals section press OK to scroll down options</i> <i>Press SELECT to select activation/deactivation of components</i> ➤ PUMP ➤ CHARGE SV ➤ DRAIN PUMP ➤ ULTRASONIC (warning: never start ultrasounds with empty tub) ➤ RS TUB (warning: never start the resistance with empty tub) ➤ RS DRYING (warning: never start the resistance without starting the FAN first) ➤ FAN ➤ EXIT
SELECT	OK	TEST CYCLE Quick test cycle
SELECT	OK	COUNTER This section shows the total cycles performed by the tool
SELECT	OK	EXIT The tool goes back to the operating mode.

		ULTRASONIC Select the operating frequency of the transducers (from 37 kHz to 50 kHz with step 1). <u>The default value is 39 kHz.</u>
		SWEEP Select the frequency difference of the ultrasonic wave during the function SWEEP (from ± 1 kHz to ± 10 kHz with step 1 including the selection NO). <u>The default value is ± 1 kHz.</u>
		SWEEP TIME Select the time in which the function of the scale wave remains flat. (operating values selectable are: 1ms / 5ms / 13ms / 25ms / 50ms / 100ms / 200ms / 500ms / 1000ms / 2000ms – more the value is low and more the frequency change will be fast). <u>The default value is 13ms.</u>
		DRY Set the duration time of the drying stage (operating values available are: 600s. / 900s. / 1200s. / 1500s. / 1800s.) <u>The default value is 900s.</u>

06

ERROR

THE APPLIANCE DOES NOT SWITCH ON	1	Power supply system check
	2	Power cord connections check
	3	Fuses integrity check
DROP OF ULTRASOUND EFFICIENCY	1	An excessive load causes a poorer efficiency of the ultrasound cleaning
PRESSING  THE CYCLE DOES NOT START	1	Wait for the end of drain
	2	Check connection of the membrane keypad to the board
	3	Replace the membrane keypad

WARNING		OPEN DOOR If the door is open while the device is operating, it automatically stops the running process.	<i>Close the door</i>
		MINIMUM TANK LEVEL DISINFECTANT / CLEANER Lack of liquid in the tank	<i>Fill the tank</i>
		MAXIMUM TANK LEVEL DISINFECTANT / CLEANER	<i>Do not add liquid</i>
		VESSEL DRAIN	<i>Wait draining stops</i>

RESET

In case of error reset the machine, hold down  for 8 seconds.

START PUMP

Start the pump for 30 second, hold down  for 8 seconds.

ERR 005

PROBLEM DURING WATER LOADING PHASE

It occurs if there are the following conditions:

If it takes more than 4 minutes

HOW TO OPERATE

Check that the water charge tap of your water system is open

Check the system pressure (1.5 bar minimum - 2 bar maximum)

Check that there are no obstructions in the inlet water pipe

Check that there is no water in the tank and then start a new cycle

Check the cleanliness and integrity of the level probe

Clean vessel filter, water inlet filter on EV charge

Check operation EV charge (enter the technical menu and enable the functioning of the single device) if necessary replaced EV charge

Check continuity in wiring that connects EV charge with Main Board and if necessary replace it

Check the WATERBLOCK (OPTIONAL) is not blocked

Replace the board (L)

ERR 007

PROBLEM DURING WATER HEATING PHASE

It occurs if the PT100 probe does not perceive an increase of 1°C in 5 minutes



HOW TO OPERATE

Replace the PT100 (K)

Replace the board (L).

Check the resistance of the silicon heaters (I) (correct medium value 150 Ω)

ERR 008		PROBLEM DURING WATER DRAIN PHASE
It occurs if there are the following conditions: <i>If at the end of the drain the tank signals MAX level</i>		
 HOW TO OPERATE <i>Remove the drain cap from the drain coupler</i> Check that there are no breaking or obstructions in the drain pipe. Replace pipe or remove obstructions Clean the filter inner the tank. Check electrical continuity and integrity of the cable connecting EV drain and Main Board. If necessary, replace the cable. <i>Replace the board.</i>		
ERR 009		PROBLEM DURING DRYING PHASE
<i>It occurs if the PT100 probe does not perceive an increase of 1°C in 1 minute</i>		
 HOW TO OPERATE <i>Replace the PT100</i> <i>Replace the board</i> Check operation of the fan and if necessary replace it Check operation of the armored resistance and if necessary replace it		

ERR 010		READING ERROR OF WATER LEVEL
During tools washing the level of probe doesn't marks MAX		
 HOW TO OPERATE		
Check the continuity in wiring and the operation of the level sensor		
Check that there is no leak in the pipes inside the machine		
Check that there is no vacuum in the drain tube		
Replace the Main Board		

ERR 011		PROBLEM DURING WATER INFLOW PHASE
It occurs if there are the following conditions: • If it takes more than 4 minutes (second water inflow)		
 HOW TO OPERATE		
Check that the water charge tap of your water system is open		
Check the system pressure (1.5 bar minimum - 2 bar maximum)		
Check that there are no obstructions in the inlet water pipe		
Check that there is no water in the tank and then start a new cycle		
Check the cleanliness and integrity of the level probe		
Clean vessel filter, water inlet filter on EV charge		
Check operation EV charge (enter the technical menu and enable the functioning of the single device) if necessary replaced EV charge		
Check continuity in wiring that connects EV charge with Main Board and if necessary replace it		
Check the WATERBLOCK (OPTIONAL) is not blocked		
Replace the board		

ERR 012		PROBLEM DURING WATER LOADING PHASE
It occurs if there are the following conditions: • <i>If the Overflow detector indicates MAX</i>		
 HOW TO OPERATE <i>Check the wiring of the sensor level (enter the technical menu and Check the sensor level, short-circuit the tank and the probe and verify that MAX signals))</i> <i>Clean or replace the sensor level</i> <i>Replace the board</i>		


CAUTION: RISK OF ELECTROCUTION

Unplug the tool from power supply before any intervention.


P1

Opening of the tool front panel

The front panel houses the membrane keypad and is internally connected by a flat cable to the control electronic board. The front panel also houses the display also connected by a multi-wire cable to the control electronic board. Be careful not to damage cables during the handling.

- 1) Unscrew the two screws located on the lower side of the panel.
- 2) Lower the panel and unscrew it.
- 3) Disconnect the flat cable connecting the membrane keypad and the control board housed inside the tool
- 4) Disconnect the multi-wire cable connecting the display to the control board housed inside the tool


P2	Opening of the plastic cover of the lid
	The plastic cover of the lid is fixed to its steel support by means of 10 screws. 1) Open the cover 2) Unscrew the 10 fixing screws 3) Close the door 4) Lift the plastic panel

P3


Removal of the rear panel

The rear panel houses all couplers for the appliance utilities and the connecting socket to the power mains. Be careful not to damage hoses and electrical cables during the handling.

- 1) Disconnect the charge and drain utilities of the tool.
- 2) Follow procedures **A1** and **A2** (removal of the front panel and of the upper cover).
- 3) Remove the cartridge filter to eliminate any water deposit from the circuit.
- 4) Disassemble from the upper fixation the gas spring after removing 3 of the 4 locking seeger rings and then extracting the fixing pin.


- 5) Unscrew the 17 fixing screws.


- 6) Disconnect the internal silicone hoses of couplers and electrical cables of the power socket, as well as the corrugated hose previously unloosened, from its seat.

P4**Removal of the power board**

1. Follow procedure **P1** to remove the front panel
2. Unscrew the 2 self-locking nuts to remove the support plate from basement

3. Then remove connectors of the board and unscrew the fixing screws.


P5

Tub removal

ATTENTION

Before removing the tub use a capacimeter to measure (at the ends of the connector on the black and white cable) the global capacity of transducers both between and 12 η F 16 η F to evaluate the actual efficiency of the transducers.

1. Follow procedure **P1** to remove the front panel
2. Follow procedure **P2** to remove the board
3. To remove the tub strongly press from the bottom to the top to unglue the tub edges from the adhesive tape fixing them.


P6

Cleaning of the filter of the charge solenoid valve

1. Close the feeding cock of the system.
2. Unscrew the ring nut of the charge hose from the rear coupler (**M**).
3. Clean the mesh filter.
4. Remount the charge hose.
5. Re-open the feeding cock.

P7**Replacement of the drain PUMP**

1. Remove the front panel and the power board following procedures **P1** and **P4**.
2. Then extract the fittings from the pump, extract fastons from contacts and remove the self-locking nuts.

3. Replace the drain pump and remount it by fixing it on the bracket at first and then carrying out the electrical and hydraulic connections.


P8

Disassembly of the charge SV

NOTE

*Before replacing the exhaust EV control the ohmic value.
Set the meter value 20k -> Ω
The correct value is 3.95 on the display*



1. Remove the front panel and the power board following procedures **P1** and **P4**.
2. Detach electric connections and the Teflon hose from the solenoid valve (the hose is attached through a John Guest quick coupler).
3. Dismount the solenoid valve by intervening on the internal nut and on the external screw.

4. Remount the new solenoid valve by following the above procedure in the opposite sequence.


P9**Sensor level replacement**

1. Remove the rear panel following procedure **P1**.
2. Remove the clamping nuts and electric.
3. Remount the level sensor new


P12**Disassembly of the peristaltic pump**

1. Remove the rear panel as indicated in procedure P1.
2. Remove the pump **open the fixing tabs as shown in the figure**, unscrew the ring nuts to extract hoses from their housing and detach electric connections.
3. Remount the new pump by following the above procedure in the opposite sequence.

REPLACEMENT OF PUMP PIPE

Remove power to the pump, so that the engine can't start accidentally.

DISASSEMBLY

- Remove the transparent protection cover through the 2 screws.
- Position the impeller with the 2 white rollers in a vertical position.
- Remove from its seat the ring on the left and then pull the pipe upward.
- Manually rotate the rotor clockwise in order to pull up even the connection put at the right of the pump.

ASSEMBLY

- Position the impeller with the 2 white rollers in a horizontal position.
- Press the pipe into his place, following gradually the path.
- Remove from its seat the ring on the left and then pull the pipe upward.
- Manually rotate the rotor clockwise in order to enter also the right connection of the pump.
- Enter again the transparent protection cover, pressing it to its place and tightening the 2 screws.

8

TESTING PROCEDURE

At the end of any assistance service is a good rule:

- Carry out a visual check to verify the correct installation or any possible leak.
- Run at least a loop to test the success of the operation.
- If necessary execute an electrical test (EN 61010)

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