AuraX Blue-22

Art. AB001ZPB

ISTRUZIONI PER L'USO	***ITALIANO***
USER INSTRUCTIONS	***ENGLISH***
MANUEL D'INSTRUCTIONS	***FRANÇAIS***
GEBRAUCHSANWEISUNG	***DEUTSCH***
INSTRUCCIONES DE USO	***ESPAÑOL***

0ZPBI0001 rev_8

22/05/2025

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CONTENTS

1. INTRODUCTION

2. USE AND INTENDED USE OF THE AUTOCLAVE

3. SYMBOLS

4. SAFETY

4.1 SAFETY MARKING

4.2 SAFETY DEVICES

4.3 SAFETY NOTES

4.4 DISPOSAL

5. TECHNICAL DATA

6. LIST OF OUTPUTS AND INDICATORS

7. UNPACKING

8. ACCESSORIES

9. INSTALLATION

10. USER INSTRUCTIONS

10.1 TURNING ON THE AUTOCLAVE AND BAROMETRIC ALIGNMENT

10.2 CONNECTIONS

10.3 DIRECT MANUAL LOADING INTO THE CLEAN WATER TANK

10.4 LOADING OF CLEAN WATER TANK THROUGH RECIRCULATION

10.5 LOADING OF CLEAN WATER TANK VIA OSMOSIS SYSTEM

10.6 WATER CHARACTERISTICS TO BE USED

10.7 LOADING MATERIALS INTO THE AUTOCLAVE

10.8 START OF STERILISATION CYCLE

10.8.1 OPERATING CYCLES

10.8.2 TEST CYCLES

10.9 CHOICE OF STERILISATION CYCLE

10.10 AUTOCLAVE ACTIVATION

10.11 START OF PROGRAMMING (START OF TIMER)

10.12 EXTRA DRYING

10.13 RUNNING THE CYCLE

10.14 END OF STERILISATION CYCLE

10.15 UNLOADING STERILISED MATERIALS

10.16 USED WATER DISCHARGE

10.17 STERILISATION CYCLE INTERRUPTION

10.18 CYCLE DIAGRAM

10.19 CYCLING DISABLEMENTS

11. DATA MANAGEMENT

11.1 VIEWING A CYCLE REPORT

12. SETTINGS

12.1 LANGUAGE

12.2 DATE AND TIME

12.3 MAINTENANCE

12.4 OSMOSIS

12.5 VALIDATION

12.6 USER DATA

12.7 RECIRCULATION

12.8 REMINDER

12.9 UNIT OF MEASUREMENT

12.10 PRE-HEATING

12.11 INFO SYSTEM

12.12 DISPLAY

13. SERVICE

14. INSTALLATION OF OSMOSIS SYSTEM

15. MAINTENANCE

15.1 ROUTINE MAINTENANCE FOR USERS

15.2 ROUTINE MAINTENANCE FOR AUTHORISED TECHNICIANS

16. SOLUTION TO OPERATIONAL PROBLEMS

16.1 THE AUTOCLAVE DOES NOT DRY PROPERLY

16.2 THE AUTOCLAVE CHAMBER HAS TURNED WHITE

16.3 THE AUTOCLAVE CHAMBER HAS GREEN-BLUE STAINS

16.4 THE STERILISER DOES NOT TURN ON

16.5 THE AUTOCLAVE DOES NOT START THE STERILISATION CYCLE

16.6 STAINS ON THE INSTRUMENTS

16.7 PRESENCE OF WATER ON THE SUPPORT SURFACE OF THE STERILISER

17. WARNING AND ALARM MESSAGES

17.1 WARNING MESSAGES

17.2 ERROR MESSAGES

18. PROCEDURE FOR SERVICE AND ASSISTANCE

19. SUMMARY OF CONSUMABLE SPARE PARTS

20. CYCLES TABLE

THE MANUFACTURER RESERVES THE RIGHT TO MAKE TECHNICAL IMPROVEMENTS TO THIS MANUAL WITHOUT PRIOR NOTICE. THIS MANUAL IS THE EXCLUSIVE PROPERTY OF TECNO-GAZ S.P.A. REPRODUCTION AND TRANSFER TO THIRD PARTIES WITHOUT AUTHORISATION IS PROHIBITED BY LAW

01 INTRODUCTION

Dear Customer,

Thank you for choosing our autoclave. We will be able to reciprocate your trust with the utmost attention and with a service that fully meets your expectations.

Before using this autoclave, please read the user manual carefully and then store it in a place accessible to all sterilisation operators.

Sterilising means adopting a precise working methodology and following specific operating protocols:

DISINFECTION: mandatory phase, to protect the safety of the designated operators, to be performed with immersion in chemical liquids or thermal disinfection;

CLEANING: the most important phase that ensures the removal of all types of residues, chemical and organic. Ultrasonic tanks are the most suitable instrument;

DRYING: an essential step, which avoids corrosion of the instruments and interference with the sterilisation cycle;

PACKAGING: an essential phase in order to maintain sterility over time;

STERILISATION: final steam sterilisation stage.

The autoclave is the key point of this methodology.

We remind you that failure to perform all the various phases of the sterilisation process may affect the final result.

For installation, maintenance and assistance, contact only authorised technicians. We invite you to use and request only original spare parts.

02 USE AND INTENDED USE OF THE AUTOCLAVE


The purpose of the autoclave is to sterilise the three types of load required by the EN13060 standard, specifically:

	22-Litre Chamber
FERROUS OR SOLID MATERIALS Cavity-free and unobstructed instruments for steam penetration	*max 5.5 kg.*
POROUS BODIES Simple materials or compounds able to absorb fluids (fabrics, gowns, gauzes, dressings, etc.)	*max 2 kg.*
CABLE BODIES Materials or devices with cavities, obstructions, etc. These are divided into two types, classified according to length and diameter. Indicatively: **TYPE A:** turbines, handpieces and devices with blind or small holes **TYPE B:** cannulas, tubes or devices with considerable passages	*max 5.5 kg.*

The autoclave must be used <u>exclusively</u> for the sterilisation of instruments and materials compatible with the steam sterilisation system and in general exclusively for the uses envisaged by the manufacturer. Always ensure that the loads subjected to sterilisation are able to withstand the temperatures and pressure of the chosen cycle.



03 SYMBOLS

	WARNING: Follow the instructions to avoid damaging the device or compromising successful operation
	DANGER: Safety instructions to avoid injury
	Pursuant to Directive 2002/96/EC, the product must not be disposed of as municipal waste at the end of its operational life
	Alternating current
	ATTENTION: high voltage. Risk of electrocution!
	Fuse
	Read the user manual carefully
	Produced on -- / -- / --
	Manufacturer
	Batch Code/Batch Number
	Catalogue Number
	Serial number

	Temperature limits
CE	CE marking
Pos.	Position (number that identifies a component in the figures or diagrams).
Fig.	Figure
	Wi-Fi Connection Active
	Maximum Clean Water Level
	Minimum Clean Water Level
	Maximum Used Water Level
	Back
	Information
	Door open / Door closed

04 SAFETY

4.1 Safety marking

	*DANGEROUS VOLTAGE*	 ATTENZIONE ATTENTION ATTENTION ACHTUNG ALTA TEMPERATURA HIGH TEMPERATURES TEMPERATURES ELEVEES HOHE TEMPERATUR	*HIGH TEMPERATURE*
 ATTENZIONE TOGLIERE TENSIONE PRIMA DI RIMUOVERE IL COPERCHIO WARNING DISCONNECT THE MAINS SUPPLY BEFORE REMOVING THIS COVER	*TURN OFF THE POWER BEFORE REMOVING THE LID*		*EARTHING CONNECTION*

4.2 Safety devices

1. Safety valve calibrated +2.4 barg +/- 0.1 bar;
2. Automatic reset safety thermostat on steam generator;
3. Manual reset safety thermostat on the heating band;
4. Mains fuses - 2 x 6.3x32 mm – F 15A 250V;
5. Electronic board fuses – 4 x 5x20 mm – T 5A L 250V + 2 x 5x20 mm – F 16A 250V;
6. Power supply protection from short circuit, overload and overvoltage;
7. Door position safety microswitch;
8. Internal pressure recovery system managed by microprocessor;
9. Continuous monitoring of internal parameters.


4.3 Safety Notes

1. The manufacturer is responsible for the product placed on the market in accordance with the current legislation. **Liability lapses** when operations are performed on the device, or on part of it, by unqualified personnel or with the use of non-original spare parts.
2. The room where the autoclave is installed must not be at potential risk of explosion and/or fire.
3. The autoclave must be installed in an environment that complies with the current legislative requirements.

4.4 Disposal

Refer to annex
Cod. 0Z00H0004

05 TECHNICAL DATA

		22-Litre Chamber
GENERAL DATA	***Medical device risk class - 2017/745/EU***	IIb
	Insulation Class - DIN EN 61140 and PN-EN 61140:2005	Class I
	Use Environment	Internal Use
	Degree of protection (IP) – IEC EN 60529 and IEC EN 60670-1	IP20
MECHANICAL DATA	***Work temperature***	+5°C ÷ +40°C
	Max altitude	2,000 m
	Relative humidity MAX at 30°C	80%
	Relative humidity MAX at 40°C	50%
	Overall dimensions *(W x H x D)*	528 x 467 x 562 mm
	Open door footprint	523 mm
	Weight *(empty tanks + empty chamber)*	approx. 60 kg .
	Weight *(full tanks + full chamber)*	approx. 74 kg.
	Maximum weight per support area (full tanks + full chamber)	17.5 N/cm^2
	Maximum sound power level	< 70 db


ELECTRICAL DATA	***Power supply voltage***	230 Vac +/-10% Single Phase
	MAX. absorbed power	2.3 kW
	Frequency	50/60 Hz
	Power cable	Cable: 3x1mm^2, 16A - 250V, 2P+E, IP20, H05VV-F3G,
		Connector 1: Schuko Angled Plug 90°C, CEE-7/VII, DIN 49441-R2
		Connector 2: Angled, IEC 60320-1/en 60320-1 Standard IEC C19 16A
	Fuses	2 x 6.3x32 mm – F 15A 250V 4 x 5x20 mm – T 5A L 250V 2 x 5x20 mm – F 16A 250V
	Transmitted heat	3.6 MJ / h
STERILISATION CHAMBER	***MAX working PRESSURE***	+2.3 barg
	MIN working PRESSURE	-0.86 barg
	MAX safety pressure	+2.4 barg
	MAX temperature	137°C
	Material	AISI 304 Stainless Steel
	Dimensions	Ø 260 x 433 mm
	Sterilisation Chamber Total Volume	approx. 22 l
	Sterilisation Chamber Useful Volume (with tray-holder inserted)	approx. 12 l

CLEAN WATER TANK	***Volume***	approx. 4 l (Water at MAX level) approx. 2.2 l (Water at MIN level)
	Cycles executable with a load (without recirculation)	from 3 to 7
	Material	PP (Polypropylene)
WATER TANK USED	***Volume***	4 l (Water at MAX level)
	Executable cycles (without recirculation)	from 5 to 16
	Material	PP (Polypropylene)
	Max. discharge water temperature	90°C
BACTERIOLOGICAL FILTER	***Diameter***	54 mm
	Filter capacity	0.027 μm
	Filter Thread	G1/8
STEAM GENERATOR	***Power***	2.2 kW
	Material	ADC12 aluminium
HEATING BAND	***Power***	1.3 kW
	Insulating material	Mica


06 LIST OF OUTPUTS AND INDICATORS

LIST OF OUTPUTS AND INDICATORS – FIG. A		
	01	Bacteriological Filter
	02	Motorised Closing System
	03	Used Water Discharge Tap (Front)
	04	Clean Water Discharge Tap (Front)
	05	Discharge Filter
	06	USB Port
	07	Main Switch
	08	LCD Screen/Control Panel
	09	Safety valve
	10	Used Water Discharge Tap (Rear)
	11	Clean Water Discharge Tap (Rear)
	12	Clean Water Overflow
	13	Power Cable Connection
	14	Demineraliser Connection
	15	Fuse


FRONT VIEW

SIDE VIEW

REAR VIEW

07 UNPACKING

The autoclave is shipped in packaging suitable for the transportation, handling and protection of the autoclave itself.

The packaging must not be knocked and must be handled with care, avoiding rolling or dropping it.

If no handling equipment is available, always ensure that two persons handle the packed autoclave.

The autoclave is supported with a wooden pallet, enclosed in corrugated cardboard and internally reinforced with cardboard compounds.

To unpack the autoclave, open the corrugated cardboard, remove the reinforcing parts and extract it using the straps supplied.


	***Handling with the use of straps must only take place performed by two persons.***
	***Never lift the autoclave from the bottom of the door. Such incorrect operation could damage the equipment.***


<u>ATTENTION:</u> Always keep the original packaging.


Inside the packaging you will find the following documents:

- *DECLARATION OF CONFORMITY:* to be kept.
- *INSTALLATION REPORT - TESTING AND WARRANTY CONDITIONS:* which you must complete when installing the machine, following the instructions indicated on the form.
- *QUICK USE GUIDE:* to be stored near the machine.
- *SAFETY VALVE DECLARATION OF CONFORMITY.*
- *ANNEXED WEEE SHEET: 0Z00H0004.*

08 ACCESSORIES



TRAY-HOLDER	
	22L-Chamber
Code	**1ZPBA0098**
Material	Stainless Steel
Dimensions *(W x H x D)*	234 x 200 x 376
Image	 ***Fig.1***
Envisaged equipment	1

NARROW TRAY	
	22L-Chamber
Code	**1ZPBA0061**
Material	Anodised Aluminium
Dimensions *(W x H x D)*	185 x 17 x 359
Image	 ***Fig.2***
Envisaged equipment	2

LONG TRAY	
	22L-Chamber
Code	**1IRAA0072**
Material	Anodised Aluminium
Dimensions *(W x H x D)*	221 x 17 x 359
Image	 ***Fig.3***
Envisaged equipment	3


TRAY EXTRACTION KEY

Use to extract and manoeuvre the trays *(chapter 10.3).*

Code	DANA008	
Image	 Fig.4	 Fig.5
Envisaged equipment	1	


CHAMBER/DOOR SEAL CLEANING SPONGE

Use it to clean the sterilisation chamber and the door seal *(chapter 16).*

Code	CPMG004
Image	 Fig.6
Envisaged equipment	1


FRONT AND REAR TAPS WATER DISCHARGE TUBE

Use to discharge water from the front taps of the autoclave (*chapter 6 - Pos. 4 and Pos.3*).
The same tube can be used to discharge the rear taps through the fitting shown in Figure 8.

Code	DANA130
Image	 Fig.7
Envisaged equipment	1


REAR TAPS DISCHARGE FITTING

Screw on the right rear tap (*chapter 6 - Pos.11*) to empty the loading tank; screw on the left rear tap (*chapter 6 - Pos.10*) to empty the discharge tank.

Code	**CPRG096**
Image	 ***Fig.8***
Envisaged equipment	1


PIPE FOR CLEAN WATER OVERFLOW DISCHARGE

Connect one end of the tube to the clean water rear overflow (*chapter 6 - Pos.12*), the other to a container for water recovery.

Code	**SXBA799**
Image	 ***Fig.9***
Envisaged equipment	1


POWER CABLE

Attach the end of the connector (*Pos. B*) to the rear panel (*chapter 6 - Pos. 13*) then the plug (*Pos. A*) directly to the power outlet of the electrical system.

Code	**3ELEE0008**
Image	 ***Fig.10***
Envisaged equipment	1


USB STICK

Use to save the autoclave cycles (*chapter 6 – Pos.6*).

Code	**3ELND0020**
Image	 ***Fig.11***
Envisaged equipment	1

M3 ALLEN KEY FOR DOOR ADJUSTMENT


Use the key to adjust the door *(chapter 16)*.

Code	**3MECL0003**
Image	 ***Fig.12***
Envisaged equipment	1

IMPORTANT
Request and use only and exclusively original accessories.


09 INSTALLATION

STERILISER SIZES

A	99 mm
B	418 mm
C	45 mm
D	467 mm
E	528 mm
F	562 mm
G	400 mm
H	1011 mm
I	122 mm
L	588 mm

COMPARTMENT DIMENSIONS

M	Minimum 767 mm
N	Minimum 628 mm
Depth	Minimum 535 mm
O	Minimum 50 mm
P	Minimum 50 mm

1. Install the autoclave in environments suitable for sterilisation.
2. The room must be adequately lit and ventilated, as required by the directives in force.
3. Install the autoclave away from heat sources and splashes of water.

4. Place the autoclave on a perfectly horizontal surface, adequately sized to support it (max.80 kg).

5. Position the autoclave at such a height that the user can inspect the entire sterilisation chamber and tanks and in order to be able to clean them easily.

6. Never lean against the door.

7. **Leave a clearance of at least 10 mm at the rear of the appliance** in order to ensure sufficient space for the electrical and hydraulic connections and for the necessary ventilation.

8. **Leave a clearance of at least 50 mm in the left side of the appliance** in order to ensure sufficient space to open the door and for the necessary ventilation.

9. **Leave a clearance of at least 50 mm on the right side of the appliance** in order to ensure sufficient space for the machine to be switched on and for the necessary ventilation.

10. **Leave a clearance of at least 300 mm at the top of the appliance** in order to be able to inspect and clean the water tanks.

11. Connect the tube supplied to discharge the clean water overflow discharge at the rear (*chapter 6 – Pos.12).*

12. Always make sure that the electrical system to which the autoclave is connected complies with the current regulations and is sized appropriately for the characteristics of the appliance.

13. Take the supplied power cable and plug the female socket into the plug of the rear panel of the autoclave (*chapter 6 – Pos.13*).

14. Connect the electrical plug to the system, making sure that it is adequate for the power supply of the machine.

15. Turn on the autoclave by pressing the main switch (*chapter 6 – Pos.7*).

16. After registering the first user, press the "*OPEN*" key and remove all the bags used to pack the individual accessories contained within the sterilisation chamber.

17. Leave only the tray-holder with the trays inside the sterilisation chamber. Place all the other accessories in an external compartment available to the operators.

NOTE:	***Avoid connecting with extensions, reductions or adapters; otherwise, micro power interruptions could be created with consequent alarm signal.***


10 USER INSTRUCTIONS

After installing the autoclave, proceed with preparation and use.

10.1 Autoclave ignition and barometric alignment


Press the main switch (*chapter 6 - Pos.7*).
Never turn on the steriliser with USB stick inserted.

Create new administrator:


Click the New User Name box and enter your text.
Click the Pin box and enter a number of at least 4 digits:


Upon opening, the autoclave will switch to the operating configuration.

10.2 Connections


Pressing the "Connections" button will display the following screen, where the autoclave's SSID is shown, along with the button to start scanning for available networks and the buttons to enable and disable the Wi-Fi connection.

On the next screen, all available Wi-Fi networks are displayed as a list. Select the network you want to connect the autoclave to and press the confirm button.

Once the desired network is selected, you need to enter the access password and press the confirm button. By pressing the "shift" icon, you can access special characters.

Successful connection will be indicated by the appearance of the Wi-Fi symbol.


If a software update is available, the above-mentioned symbol will light up green. To perform the software update procedure, go to the INFO SYSTEM section of the settings menu (see Chapter 12.11).

THE AUTOCLAVE IS NOW READY FOR USE.

If the machine is started with the door open, barometric alignment will occur automatically, and the autoclave will proceed directly to the operational configuration.

10.3 Direct manual loading into the clean water tank

Lift the cap on the top of the autoclave, on the left side.

With the machine switched on, introduce water into the clean water tank until the "max" symbol with yellow LEDs is switched on.

10.4 Loading of clean water tank through recirculation

After manually loading the tank, activate the water recirculation (*chapter 13.7*); in this way, whenever the maximum level in the dirty water tank is reached, a recirculation pump will automatically be activated which, by means of the appropriate filters, consisting of resins and activated carbons, will return the distilled water to the clean water tank.
Thanks to this system, the user will not need to empty the used water tank every time or refill the clean water tank.

10.5 Loading of clean water tank via osmosis system

Connect the Pura Plus system as indicated in the relevant user manual and activate the osmosis system (*chapter 13.4 and chapter 15*).


10.6 Water characteristics to be used


TABLE OF QUALITY LEVELS ESTABLISHED BY THE EN13060 STANDARD			
Evaporation residue	≤	**10**	mg/l
Silicon oxide	≤	**1**	mg/l
Iron	≤	**0.2**	mg/l
Cadmium	≤	**0.005**	mg/l
Lead	≤	**0.05**	mg/l
Remains of heavy metals except iron, cadmium, lead	≤	**0.1**	mg/l
Chloride (Cl')	≤	**2**	mg/l
Phosphate (P20s)	≤	**0.5**	mg/l
Conductivity (at 20°C)	≤	**15**	µs/cm
Ph value (degree of acidity)	**5 ÷ 7.5**		
Appearance	Colourless, clear, no deposits		
Hardness (and alkaline earth ions)	≤	**0.02**	mmol/l


10.7 Loading materials into the autoclave

Arrange the materials to be sterilised on the trays supplied, being sure to:

- arrange the objects sufficiently spaced apart from each other;
- never overlap or stack materials;
- arrange the packed instruments, always with the paper part facing upwards;
- never place the materials in contact with the sterilisation chamber or with the closing door;
- arrange pincers and scissors with the blades open;
- arrange the tray-holder in a perfectly horizontal position.

The tray holder for surgical cassettes can be used by rotating it 90° from its normal position. When the procedure is completed, the display will show the word “OPEN”.

10.8 Start of sterilisation cycle

After the steps indicated in chapter 10.6, choose the sterilisation program most suitable for the load prepared.

The autoclave is equipped with two series of cycles, represented in different colours:

A - Operating cycles
B - Test cycles

10.8.1 Operating cycles

A - All the operating cycles involve a fractional vacuum system (3 vacuums) and can sterilise hollow, porous and solid bodies; however, they differ according to the priority of the user and the characteristics of the materials to be sterilised. The possible sterilisation temperatures are 121°C and 134°C.

- Universal 134°C cycle: it is used for all materials that are not thermolabile up to a load of 5.5 Kg if they are solid or up to a load of 2 Kg if they are porous; the cycle is suitable for wrapped or non-wrapped loads.

- 121°C Universal cycle: it is used for thermolabile or otherwise sensitive materials, up to a load of 5.5 Kg if they are solid or up to a load of 2 Kg if they are porous; the cycle is suitable for wrapped or non-wrapped loads.

- 134°C Flash cycle, (**Class S autoclave cycle)** created to meet the immediate needs of the instrumentation, it has a short cycle execution time due to the exclusion of the drying phase; it maintains the characteristics of the 134°C Universal cycle. Not suitable for wrapped loads.

- 134°C Prion cycle: designed for Creutzfeldt-Jakob disease (mad cow disease), it has a load capacity of 5.5 Kg if the material is solid or up to a load of 2 Kg if the material is porous; the cycle is suitable for wrapped or unwrapped loads. The exposure time is much longer than that of a 134°C Universal cycle.

Please refer to the table in Chapter 21 for further details.

10.8.2 Test Cycles

B - The available test cycles are:


- Vacuum Test;
- Helix/B&D Test

- Vacuum Test: used to verify the correct seal of the sterilisation chamber by measuring the change in depression over a period of 15 minutes.
This cycle must be performed with the sterilisation chamber empty, with only trays and tray-holders.
The cycle is not started if the temperature in the chamber is too high as the pressure is affected by high temperatures. For this reason it is advisable to perform the test at the beginning of each working day when the chamber is at room temperature.

- Helix/B&D Test: this test is identical to a 134°C Universal cycle, but with a sterilisation duration of 3.5 minutes instead of 4.
By means of an appropriate instrument (Test Helix) it is possible to verify the penetration of steam inside hollow bodies.
This test is also useful to evaluate the correct penetration of steam inside porous bodies (Bowie&Dick Test).

10.9 Choice of sterilisation cycle

After choosing the program, the following screen will appear:

The following information will appear in the top right:

- Total cycle time;
- Sterilisation time;
- Nominal pressure in sterilisation;
- Maximum Solid Load;
- Maximum Porous Load

ATTENTION:

The drop-down menu on the left, after 10 seconds, will compact to make room for the main screen. Press on the left side of the screen on any of the icons to bring up the various menu items again.

10.10 Autoclave activation

If the autoclave is not activated, the cycle start is inhibited and the activation screen appears.
Connect to the website https://s.tecnogaz.com.
Go to the “End Users” section and authenticate.
In the “Your Devices” section, under the stereotyped image of the device, click the key:

The “AURA ACTIVATION CODE” panel opens:


Write the ID code found on the autoclave and click on "generate key" the relative activation code will appear to be inserted in the appropriate box.

Proceed with activation by pressing the green button.

10.11 Start of Programming (Start of Timer)


Each cycle and selectable test provides the possibility of programming the start time through the "Start Timer" icon.

Select the parameter to be modified (hour, minute, day, etc.) and modify the value with the arrows. Once the cycle start date has been set, press the lock icon that will become closed when opened.
At this point, press the "Back" key.

In the main screen of the cycle the alarm will have turned yellow, to indicate the fact that a Start Timer has been set.


Press the "Start" key; the screen below will appear where the user can decide whether to stop the cycle programming through the "stop" key or instantly start the program.
Leaving the screen active, the cycle will automatically start at the set time and date.

NOTE:	***Before setting the start of the programmed cycle, be sure to have closed the door correctly.***

10.12 Extra Drying

The 134°C Universal, 121°C Universal and 134°C Prion cycles provide the possibility of increasing the drying period through the "Extra drying" key. This option is used to extend the phase period from 0 to 15 minutes.

Press the arrows to lengthen or shorten the extra drying time. Once the extra drying time has been set, press the lock icon that will become closed when opened.
After determining the desired minutes, press the "Back" key.


In the main screen of the cycle, the icon relating to the extra drying will have turned blue, to indicate the fact that an extra drying has been set in the cycle to come.

After pressing the "start" key, a cycle will begin with the duration of the drying phase increased by the established value.


10.13 Running the cycle

During the cycle, the display will show all the main parameters relating to the cycle performed; in particular, the following will be displayed:

- Progress Phase;
- Temperature in the sterilisation chamber;
- Pressure in the sterilisation chamber;
- Residual time at the end of the cycle (from the sterilisation phase onwards).


Through the "Information" button (not only during a cycle, but at any time when the machine is switched on) it is possible to access the list of additional operating parameters such as: Theoretical Temperature, Steam Generator Temperature, Heating Band Temperature, Water Used, Water Conductivity, etc.


10.13 End of sterilisation cycle


The cycle completed message will appear on the display with green LEDs attached.

If cycle validation is active (chapter 12.5), a screen will appear where you can select the user who performed the sterilization cycle along with their respective pin.

The registered user will proceed to validate the cycle by pressing the green button if they consider the load to be compliant.

10.15 Unloading sterilised materials

Wear suitable personal protection equipment in accordance with the current regulations on safety and hygiene at work. Extract the trays using the appropriate key provided (*chapter 8 - Fig.4/5)*, leave the instruments to condition and store them in environments where they cannot be contaminated.

10.16 Used water discharge

When the "Maximum Used Water Level" symbol appears, if the recirculation is deactivated (*chapter 13.7)*, the used water tank must be emptied.

Unless otherwise instructed, the operation of the autoclave is inhibited.

Take the supplied tube (*chapter 8 - Fig.7*) and insert it into the used water discharge fitting located on the front of the autoclave (*chapter 6 – Pos.3*). Place the other end of the tube in a container of at least 4 litres and unscrew the ring nut anti-clockwise; the water will fall into the container, emptying the tank.

IMPORTANT:

A - The tube housed in the collection container must never touch or be immersed in the discharged water.


B - Always wait until the discharge water has completely drained out. The maximum water level used symbol turns off when there is still water in the tank, so do not use it as a reference for this operation.


C - Any residues inside the discharge tank are due to the normal operation of the autoclave and do not create problems for the sterilisation cycles.

Once discharging is completed, tighten the ring nut and remove the tube.

10.17 Sterilisation cycle interruption

A sterilisation cycle can be intentionally interrupted by pressing the "STOP" button. The steriliser will immediately start an internal pressure reset cycle. At the end of this process, the alarm message A100 will appear on the display.

To reset the alarm, press and hold the "reset" key for at least 5 seconds.
After a manual interruption of the program, the load cannot be considered sterilised.

10.18 Cycle diagram

The image below shows the typical pressure trend inside the chamber during a cycle with standard drying.


	***CYCLE DIAGRAM*** ***T1, T2, T3 Pre-heating and vacuum fractionation phases*** ***T4 Sterilisation*** ***T5 Drying***


Fig. 21

TABLE OF TYPE TESTS AS PER STANDARD EN13060				
Type test	**134°C Universal**	**121°C Universal**	**134°C Flash**	**134°C Prion**
Pressure dynamics in the sterilisation chamber	**X**	**X**	**X**	**X**
Air leakage	**X**	**X**	**X**	**X**
Vacuum chamber	**X**	**X**	**X**	**X**
Solid load	**X**	**X**	**X**	**X**
Small porous load	**X**	**X**		**X**
Full porous load	**X**	**X**		**X**
Cable load B	**X**	**X**	**X**	**X**
Cable load A	**X**	**X**	**X**	**X**
Multiple packaging	**X**	**X**	**X**	**X**
Solid load drying	**X**	**X**		**X**
Porous load drying	**X**	**X**		**X**

TECNO-GAZ is at your disposal to provide you with sterilisation tests for the management of your sterilisation routine


10.19 Cycling Disablements

By contacting an authorized technician, it is possible to deactivate some operating cycles to prevent them from starting. The icons of the deactivated cycles are clearer like the one shown in the figure below and are not clickable.


11 DATA MANAGEMENT

The "Data Management" menu is used to view, print, export and/or delete the cycles performed by the steriliser.


By pressing the print key, one or more selected cycles can be printed.
To change the expiry date and the number of labels to be printed according to requirements, select the parameter to be modified (day, month, year, number of copies) and change its value with the arrows.

By inserting a USB stick in the appropriate port (*chapter 6 – Pos.6)* and pressing the export key, one or more cycles can be transferred to USB.
The USB stick must be formatted in FAT32.
After selecting the cycles to be exported, and after pressing the "Export" key, a confirmation message will appear; approve and wait for the data transfer.
By pressing the delete key, one or more cycles can be deleted.
Before printing, before exporting the data to a USB stick or **before deleting it**, select from the screen the cycles to be exported/deleted. In order to facilitate this operation, the user can select through the top row:

- The last 10 cycles;
- The last 50 cycles;
- All cycles.

In particular, the display screen of the cycles shows:

- Cycle number;
- Cycle type;
- Date and time of the end of the cycle;
- Cycle outcome:
 - A. Validated cycle;
 - B. Cycle interrupted.

WARNING: You can select and export a maximum of 20 cycles at a time.

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- **Start**: Moment of cycle start, Beginning of Phase 1
- **SV1 (Start Vacuum 1)**: Moment when Vacuum 1 starts, Beginning of Phase 2
- **EV1 (End Vacuum 1)**: Moment when Vacuum 1 ends, Beginning of Phase 3
- **EP1 (End Preheating 1)**: Moment when Preheating 1 ends, Beginning of Phase 4
- **EV2 (End Vacuum 2)**: Moment when Vacuum 2 ends, Beginning of Phase 5
- **EP2 (End Preheating 2)**: Moment when Preheating 2 ends, Beginning of Phase 6
- **EV3 (End Vacuum 3)**: Moment when Vacuum 3 ends, Beginning of Phase 7
- **EPH (End Preheating)**: Moment when Preheating 3 ends, Beginning of Phase 8
- **SS (Start Sterilization)**: Moment when Sterilization starts, Beginning of Phase 9
- **ES (End Sterilization)**: Moment when Sterilization ends, End of Phase 9
- **DR (Drain)**: Moment when the draining phase ends, End of Sub-phase 10-4
- **SD (Start Drying)**: Moment when drying starts, Beginning of Phase 11
- **ED (End Drying)**: Moment when drying ends, End of Phase 11
- **End**: End of Phase 12

12 SETTINGS

The settings menu is divided into two pages. It includes a wide range of machine configurations; in particular, the preferences that the user can customise are:


- Language;
- Date and time;
- Maintenance;
- Osmosis;
- Validation;
- User Data;
- Water Recirculation;
- Reminder;
- Unit of Measurement;
- Pre-heating;
- Info System;
- Display.

12.1 Language

Select the desired language.

12.2 Date and Time


To change date, time and UTC, select the parameter to be changed (day, month, year, hours, minutes, etc.) and change its value with the arrows. It is also possible to choose the date format. Then press the confirmation key.


12.3 Maintenance

Through this menu it is possible to reset the meters, once maintenance has been performed, of:


- Bacteriological Filter;
- Door Seal;
- Osmosis Filters;
- Recirculation Filters.

The autoclave is equipped with a warning system that alerts the operator to the need for maintenance to ensure future correct use.
Within one of the 4 sub-menus it is possible to reset the meter relating to that specific maintenance by pressing the "reset" key.


12.4 Osmosis


Select "Activate" and then the confirmation key, in the case of an osmosis system to be integrated into the autoclave.

12.5 Validation

Select "Activate" and then the confirmation key to validate or not the success of the sterilisation cycle at the end of the process. The conformity of the load will then be recorded in the sterilisation report.

12.6 User Data

Through this menu it is possible to delete or add a user, as well as change their password. The creation of a user is linked to the identification when the machine is switched on.


ATTENTION:

After 3 incorrect pin entries, during power-up, a message is generated requesting an unlock code. Enter:

123456789


Then go to the User screen and create a new user.

12.7 Recirculation

Activating this option will activate the recirculation and filtering system whenever the maximum level of the dirty water tank is reached. The pump will be operated until the maximum level of the clean water tank is reached.


12.8 Reminder

This function is used to display, based on the chosen days, a message that refers to the execution of the relative Test.


The reminder will appear at the same time it was activated or when the steriliser was switched on if the machine was switched off at that time.

12.9 Unit of measurement


By pressing the "Unit of Measurement" icon it is possible to set the desired units of measurement (temperature, water conductivity and pressure).

12.10 Pre-heating

Through the "Pre-heating" icon it is possible to set the maximum time (expressed in hours) in which the sterilisation chamber is kept warm (and therefore ready) for an imminent sterilisation cycle. The user can deactivate this function or vary the period of time in which it is active from a minimum of 1 hour to a maximum of 8 hours.
By default, a time of 4 hours is set.

12.11 Info System

In this screen you can view the sterilizer features such as:

- Serial Number;
- Model;
- FW Version (;
- LCD FW Version;
- Key-ID.

If the sterilizer is connected via Wi-Fi, and if new software updates are available, the WI-FI symbol icon and the “Software Update” button will light up green.


To proceed with the update, click the “Software Update” button.

12.12 Display

By pressing the "Display" icon it is possible to set the brightness of the Display; once set, press the confirmation key.

13 SERVICE

Access to this menu is reserved only for technicians authorised via PIN.

14 INSTALLATION OF OSMOSIS SYSTEM

Setting up the demineralisation system

The autoclave is designed to load the demineralised water through an external osmosis demineralisation system (*optional accessory*).

Before installing the system, the operator must program the autoclave following the instructions below:

Turn on the autoclave by pressing the main switch (*chapter 6 – Pos.7*).

When the autoclave has been brought to the operating screen, access the OSMOSIS menu (*chapter 13.4*) and press "Activate".

Exit the menu and return to the home screen.

NOTE	***With connection by means of demineraliser, if the maximum water level has not been reached, the operation of the autoclave will be inhibited.***

ATTENTION

When the maximum number of refills with the osmosis system is reached, a message on the display will inform the user. It is necessary to remember to reset the cycle meter when changing the filter on the osmosis system (*chapter 13.3*).

Connection of demineralisation systems

- Switch off the autoclave if it is switched on (*chapter 6 – Pos.7*);
- Close the tap located upstream of the demineralisation system and the clean water discharge tap (*chapter 6 – Pos.11*);
- Install the demineralisation as indicated in the manual of the demineraliser itself;
- Wrap the male thread of the tube fitting with Teflon or another component that guarantees water tightness;
- Screw the tube fitting onto the female thread of the clean water discharge (*chapter 6 – Pos.11*);
- Insert the tube exiting the demineraliser into the tube-holder fitting now screwed to the autoclave;
- Insert the demineraliser plug into the socket (*chapter 6 – Pos.14*) in the back of the autoclave;
- Open the tap located upstream of the demineralisation system and then the clean water discharge tap (*chapter 6 – Pos.11)*;
- Check that there are no water leaks;
- Switch on the autoclave;
- Perform one or more sterilisation cycles to check the operation of the connection made and especially to check for leaks.

	***At the end of the day, always close the tap located upstream of the demineralisation system***
	***Connect the demineralisation systems only to prepared autoclaves***
NOTE	***For the connection of demineralisation systems to autoclaves, also refer to what is indicated in the manual of demineralisation systems***
NOTE	***The first filling may take longer than is usual. Subsequently, the autoclave will fill automatically during the cycles.***

15 MAINTENANCE

15.1 Routine maintenance for users

Correct maintenance of the autoclave ensures efficient operation of the same and guaranteed savings in terms of time and costs linked to assistance. The following operations are to be considered mandatory and can be performed by operators.

ATTENTION

All the operations described below must be performed with the machine off and cold.

Cleaning the sterilisation chamber, trays, tray-holders and door seal	20 cycles or once a week

Periodically clean the chamber, trays, tray-holders and door seal, removing any deposits or debris; this will prevent materials from entering the discharge circuit that could create obstructions. For good cleaning use only demineralised water and the supplied sponge (*non-abrasive part – chapter 8 - Fig.6*).

Always to be performed with the machine off and with the chamber cold to avoid burns - Never use solvents, detergents, chemical solutions, descalers or other similar products.

The formation of white stains at the base of the inner wall of the chamber means that demineralised water of poor quality has been used.

Cleaning of external surfaces	Once a week

The use of ethyl alcohol diluted to a maximum of 50% with water is recommended for the cleaning and disinfecting of external surfaces. Apply the product with soaked cloth, then dry.

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For the cleaning and disinfection of external surfaces it is advisable to use a disposable cloth/paper soaked in up to 70% ethyl alcohol with water.

ATTENTION

Do not use 70% alcohol for the disinfection of other plastic surfaces.

Empty the tanks of clean and used water, remove any deposits on the bottom, wipe the cloth soaked in 70% ethyl alcohol and clean thoroughly.

Changing the Bacteriological Filter	200 cycles or when it changes to a dark colour


Replace the bacteriological filter (*chapter 8 - Fig.1*) by turning it anti-clockwise to unscrew it and clockwise to tighten it. *Only use original filters (replacement cod. DAVA101).* Once the operation has been performed, remember to reset the cycle meter (*chapter 13.3*).

Filter change (activated carbon filter and resins)	200 recirculations or when the water conductivity exceeds the limit after a recirculation (Error A301)


Only use original filters (replacement cod. AB800ZPB):

To replace the Filters:

1. Completely empty the clean water tank (chapter 6 – Pos. 4, chapter 6 – Pos 11);
2. Lift the upper lid;
3. Remove the pair of filters from the clip, taking care not to strike the level switches in the tank;

4. Disconnect the rubber tube from the metal fitting;
5. Connect the rubber tube to the replacement kit;
6. Hook the two filters together with the double clips, being sure to position the filters as shown in the figure below;
7. Hook the larger filter to the machine by pushing it towards the front of the machine until it touches the tank walls;


ATTENTION

Make sure that the assembled spare kit does not interfere with the movements of the level switches in the tank.

8. Refit the upper lid;
9. Manually fill the tank with clean water (*chapter 10.2*);
10. Reset the "Recirculation filter" meter (*chapter 13.3*)

Door adjustment	every 2 months

To preserve the functional integrity of the machine, the closing pressure of the door must be adjusted. Use the M3 Allen key supplied (*chapter 8 – Fig.14*) and loosen the screw shown in the figure (it is not necessary to unscrew it completely, simply unscrew it 3 turns anti-clockwise).

Turn the door clockwise by **a quarter of a turn** to reduce the closing pressure, or turn it anti-clockwise by **a quarter of a turn** to increase the closing pressure. Tighten the previously loosened screw to secure the new door position.

16.2 Routine maintenance for authorised technicians

TWICE-YEARLY CHECK OR EVERY 1000 CYCLES	Clean the chamber
	Clean the metal filter inside the chamber
	Replace the chamber door seal (1IRAA0020)
	Bacteriological filter replacement (DAVA101)
	Replacement of solid impurities filter in the clean water tank and solid impurities filter in the recirculation system (DARA054 x 2)
	Cleaning of brass Y-filter downstream of the radiator
	Clean the tanks
	Clean the metal filter in the clean water tank
	Check the tank level probes
	Clean the radiator and cooling fans
	Electro/pneumatic circuit integrity check
	Check that there are no traces of first disinfection liquid inside the pneumatic circuit.
	Safety valve check
	Vacuum pump performance check
	Validation in compliance with and at intervals defined by the local regulations
	Electrical safety tests every time the machine is operated
	Water Conductivity……. μS/cm
ADDITIONAL CHECK AT 2000 CYCLES	Replace valves and vacuum pump o-rings – Cleaning/integrity check of Knf membranes (CPGM025 x4 – CPGM043 x2) - Thomas (CPGM050 x4 – CPGM051 x8 – CPGM052 x2)
	Check the three N/C solenoid valves
	Check the heating band
	Check the steam generator
ADDITIONAL CHECK AT 3000 CYCLES	Replace the three N/C solenoid valves (CEEG021 x3)
	Replace the heating band (1ZPBA0021)
	Cleaning/checking the integrity of the vacuum pump membranes

16 SOLUTION TO OPERATIONAL PROBLEMS

In many cases, some alarms or errors are caused by lack of attention or knowledge of some technical and operational aspects. Below are some cases of anomalies with the related solutions.

16.1 The autoclave does not dry properly

- An excessive load has been placed inside the autoclave. Check that the load does not exceed the maximum permissible limits (*chapter 2 and chapter 21*);
- Non-original trays were used, or those of different material, without holes or with different holes. *It is advisable to use only original trays;*
- The instruments have not been arranged correctly (indications in *chapter 10.6*);
- Chamber discharge filter clogged (*chapter 16*).
- Increase the drying time using the "Extra-drying" function (*chapter 10.11*).

16.2 The autoclave chamber has turned white

- Change the type of water being used, use demineralised or distilled water, as specifically indicated in the previous chapters, and then clean the chamber.
- The whitish colour may be a consequence of the evaporation of organic materials, present on the instruments. Subject the instruments to a more suitable and thorough cleansing action.
- Check the demineralisation system installed, if any.

16.3 The autoclave chamber has blue-green stains

- There was no proper rinsing of the instruments after the cleaning phase. Rinse the instruments more carefully and thoroughly. If stains are evident, request technical assistance by telephone;

16.4 The steriliser does not turn on

- Check if the autoclave is connected to the mains with extensions, reductions or adapters; if necessary, remove these accessories and connect the autoclave directly to the electrical outlet;
- The plug of the power cable is not properly inserted into the socket;
- Absence of voltage at the power socket;
- The mains fuses are broken.

16.5 The autoclave does not start the sterilisation cycle

- The clean water tank is empty. There is a minimum level icon. Fill the tank with demineralised water;
- The used water tank is full. There is a maximum used water level icon. Discharge the used water.

These two blocks at the start of the cycle are implemented for safety reasons in such a way as to avoiding operating the steam generator without water as well as avoiding the risk of water leaking from the used water tank.
Both cases are reported by a warning *(chapter 18).*
The warnings can be consulted in *chapter 18.1* and in particular these two cases refer to warnings A005 and A004.

16.6 Stains on instruments

- The instruments turn yellow which means that a residue of chemical liquid has set on the instruments with the heat. Adequate rinsing was not performed.
- The sterilisation chamber has yellow stains. Instruments were placed in the chamber with the presence of chemical liquid that, when falling, became fixed thanks to the heat. Adequate rinsing was not performed.
- The instruments have whitish stains. Rinsing has been performed with very calcareous water and the instruments have not been dried. As a last rinse, it is advisable to use demineralised water and to dry the instruments thoroughly.
- The instruments have blackened. This is due to the fact that the instruments have strong carbon components inside.

16.7 Presence of water on the support surface of the steriliser

- Leakage of steam/water from the door seal. At the end of the cycle clean the seal and the closing sightglass with a clean cloth; check for any damage to the seal; if the problem persists, an adjustment of the door becomes necessary (*chapter 16*).

17 WARNING AND ALARM MESSAGES

If an abnormal condition occurs during operation of the steriliser, a warning or alarm message is generated.

- Warning messages are alerts generated by the system that can be resolved by the user; when they appear, the LED bar turns yellow.
- Alarm messages, instead, are related to a problem that can generally be restored by a specialist technician; when they appear, the LED bar turns red.

Warning and alarm messages are highlighted by an alphanumeric code, consisting of a letter and 3 digits.

17.1 Warning Messages

WARNING CODE	WARNING DESCRIPTION	SOLUTION
A003	Maximum Water Level in the Load Tank	Stop Water Loading.
A004	Maximum Water Level in the Discharge Tank	Empty the Discharge Tank.
A005	Minimum Water Level in the Load Tank	Fill the Clean Water Tank.
A006	Bacteriological Filter Replacement	Replace the bacteriological filter as it is depleted (*chapter 16*).
A007	Door Seal Replacement	Replace the door seal (*chapter 16*).
A008	Osmosis System Maintenance	Perform routine maintenance on the osmotic system (*chapter 15*).
A009	Recirculation System Maintenance	Replace the recirculation filters as they are depleted (*chapter 16*).
A010	Internal Maintenance	Contact Assistance to perform routine internal maintenance on the machine.
A011	No Water Recirculation Taking Place	Replace the recirculation filters as they are depleted (*chapter 16*).
A012	Low Water Quality	Empty the clean water tank and refill it with demineralised water.

A013	Temperature Too High for Testing	Go to the Vacuum Test menu with the door open and wait for the chamber to cool down.
A014	Vacuum Test Reminder	Start a Vacuum Test or ignore the message.
A015	Helix&Bowie Reminder	Start a Helix&Bowie Test or ignore the message.
A016	Close the Door	Approach the door and wait for the start of the cycle.
A017	Incorrect Pin	Access the next screen and enter the unlock code (*chapter 13.6*)
A018	Open the Door	Press "Open" on the display
A019	Recirculation not working	Check that the maximum sensor of the dirty water tank is not blocked
A020	Low water level at start-up	Fill the clean water tank with demineralised water

17.2 Alarm Messages

*If one of the alarms reappears after a short time, contact **technical assistance**.*

If an alarm message is displayed during a cycle, it is considered UNSUCCESSFUL: all the preparation and sterilisation operations must be repeated.

To reset the alarms, press and hold the "RESET" button for 5 seconds.

ALARM CODE	ALARM DESCRIPTION	SOLUTION
A100	Cycle Interrupted	Restart the cycle.
A110	Vacuum not reached in Vacuum Test	Check the condition of the seal; Check the adjustment of the door; Restart the cycle.
A111	Vacuum not maintained on first Vacuum Test phase	Check the condition of the seal; Check the adjustment of the door; Restart the cycle.
A112	Vacuum not maintained on second Vacuum Test phase	Check the condition of the seal; Check the adjustment of the door; Restart the cycle.
A113	High Temperature Change in Vacuum Test	Cool the machine with the door open; Restart the cycle.
A120	Initial heating error	Restart the cycle.
A121	Vacuum not reached	Check the condition of the seal; Check the quantity of load introduced; Check the adjustment of the door; Restart the cycle.
A122	Steam generator error	Restart the cycle.

A123	Steam generation error	Check the quantity of load introduced; Restart the cycle.
A124	Pre-heating error	Check the quantity of load introduced; Restart the cycle.
A125	Final pre-heating error	Check the quantity of load introduced; Restart the cycle.
A126	Sterilisation temperature above 124°C	Allow the chamber to cool; Restart the cycle.
A127	Theoretical sterilisation temperature above 124°C	Allow the chamber to cool; Restart the cycle.
A128	Sterilisation temperature above 137°C	Allow the chamber to cool; Restart the cycle.
A129	Theoretical sterilisation temperature above 137°C	Allow the chamber to cool; Restart the cycle.
A130	Sterilisation temperature below 121°C	Allow the chamber to cool; Restart the cycle.
A131	Theoretical sterilisation temperature below 121°C	Allow the chamber to cool; Restart the cycle.
A132	Sterilisation temperature below 134°C	Allow the chamber to cool; Restart the cycle.
A133	Theoretical sterilisation temperature below 134°C	Allow the chamber to cool; Restart the cycle.
A134	High probe reading difference	Allow the chamber to cool; Restart the cycle.

A135	Equilibration time too long	Allow the chamber to cool; Restart the cycle.
A136	Pressure discharge too slow	Allow the chamber to cool; Clean the chamber filter; Restart the cycle.
A137	Pressure discharge too slow after sterilisation	Allow the chamber to cool; Clean the chamber filter; Restart the cycle.
A138	High chamber pressure	Allow the chamber to cool; Restart the cycle.
A139	Drying duration exceeded time limit	Allow the chamber to cool; Restart the cycle.
A200	Closure exceeded maximum time	Switch off the machine; Restart the machine to open the door.
A201	Maximum absorbed current error	Switch off the machine; Restart the machine to open the door.
A202	Door signal errors	Switch off the machine; Restart the machine to open the door.
A300	Level Switches Error	Switch off and restart the machine.
A301	Water Quality Beyond the Operational Limits	Empty the clean water tank and refill it with demineralised water.
A900	Power supply interruption	Check the power supply system of the machine and of the room; Restart the cycle.
A901	SD Card missing or malfunctioning	Switch off and restart the machine.

A400	Pressure Sensor Short Circuited	Switch off the machine
A401	Pressure Sensor Not Connected	Switch off the machine
A402	Temperature Sensor 1 Short Circuit	Switch off the machine
A403	Temperature Sensor 1 Not Connected	Switch off the machine
A404	Temperature Sensor 2 Short Circuit	Switch off the machine
A405	Temperature Sensor 2 Not Connected	Switch off the machine
A406	Temperature Sensor 3 Short Circuit	Switch off the machine
A407	Temperature Sensor 3 Not Connected	Switch off the machine
A408	Serial Number Does Not Match	Switch off the machine
A409	All_Info_Alarm_Triac_Power	Switch off the machine
A410	All_Info_Alarm_Triac_Aux	Switch off the machine
A411	All_Info_Alarm_Memory	Switch off the machine

18 PROCEDURES FOR SERVICE AND ASSISTANCE

In case of a fault, service or validation, contact ***TECNO-GAZ S.p.A.***

Assistance will be provided to assess the need to return the product to the manufacturer or for the intervention of a technician and, having viewed the machine, to draw up an expense estimate, which will be forwarded to the distributor customer who will send it to the end customer, for acknowledgement and signing.

After receiving the quote signed for acceptance, the autoclave will be designated for repair and will be shipped respecting the relevant times, indicated on the quote form.

If the autoclave must be shipped to the manufacturer, follow the following mandatory instructions:

1. Use the original packaging. If this is no longer in your possession, use other appropriate packaging. The goods travel at the risk of the sender.

2. Only ship the autoclave (do not insert any components contained in the accessory kit).

3. Thoroughly clean the sterilisation chamber and the autoclave in its entirety before shipping. If it arrives dirty and with residues, the autoclave will be sent back unrepaired or will be subjected to cleaning and disinfection.

4. Always discharge the demineralised water tank.

5. Always discharge the used water tank.












6. Indicate in writing, and insert in the packaging a document where the anomaly found or the service you intend to benefit from is precisely indicated.

7. Always ship free port, otherwise the transport costs incurred will be charged.

All non-original packaging received will be disposed of.

The machine will be sent back with original packaging (the cost of the packaging will be charged to you) by the customer's shipping agent.

19 SUMMARY OF CONSUMABLE SPARE PARTS

	***TRAY-HOLDER (AuraX Blue 22 only)***		***NARROW TRAY (AuraX Blue 22 only)***
	1ZPBA0098		*1ZPBA0061*
	***LONG TRAY (AuraX Blue 22 only)***		***TRAY EXTRACTION KEY***
	1IRAA0072		*DANA008*
	***BOILER FILTER***		***SPONGE***
	DXBA091		*CPMG004*
	***REAR TAPS FITTING***		***DOOR SEAL***
	CPRG096		*1IRAA0020*
	***BACTERIA FILTER LOGICO***		***TANK FRONT DISCHARGE TUBE***
	DAVA101		*DANA130*
	***POWER CABLE***		***ALLEN KEY M3 DOOR ADJUSTMENT***
	3ELEE0008		*3MECL0003*
	***CLEAN WATER OVERFLOW DISCHARGE TUBE***		***RECIRCULATION FILTER UNIT REPLACEMENT KIT***
	SXBA799		*AB800ZPB*
	***USB STICK***		
	3ELND0020		

CHEMICAL SUPPLEMENT	*200/S*	**BOWIE & DICK TEST**	*TS001BED*
BIOLOGICAL INDICATOR	*TS002ZBK*	**HELIX TEST**	*TS004ZHT*

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