







It is recommended to check the icons on the screen, the expiration date of the electrodes and to open the case at least once a week or before use.



	Battery status	OK
	Ready to use	OK



Person authorized and trained to use the defibrillator



.....

Date	Check	yes /no	Signature
.....		
			
	Have the pads expired?		
	Open box message: Functional defibrillator		

Date	Check	yes /no	Signature
.....		
			
	Have the pads expired?		
	Open box message: Functional defibrillator		

Date	Check	yes /no	Signature
.....		
			
	Have the pads expired?		
	Open box message: Functional defibrillator		

Date	Check	yes /no	Signature
.....		
			
	Have the pads expired?		
	Open box message: Functional defibrillator		

Date	Check	yes /no	Signature
.....		
			
	Have the pads expired?		
	Open box message: Functional defibrillator		